

CALCIUM

Spectrophotometric
ARSENAZO III

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instruments: EXPRESS 550 / Plus

Reagent preparation

Reagent is ready to be used.

Instrument settings

Test Name	CA	Test	CA		
Test Bar Code					
Test Type	EndPoint	Curve Type	Blanked Linear		
Units	mg/dL	Nr. of Decimal Places	1		
Primary Wavelength	600	Secondary Wavelength			
Read Time/Interval	120	Sample Blank	No		
Factor					
Calibration Interval	168 hours	Nr. of Replicates	2		
Nr. of Calibrators	2				
Low Blank A Limit	-0.100	High Blank A Limit	2.100		
Low A Limit	-0.100	High A Limit	2.500		
Low Normal	9.0	High Normal	10.7		
Linearity Limit	18.0	Curve S.D. Limit	10.0		
Sample Volume	5 µL	Test	CA		
Rerun Dilution Ratio	1				
Sample Diluent Bottle Type	*	Predilution Ratio	1		
Reagent Diluent Bottle Type	*				
	Reagent Volume	Bar Code	Diluent Volume	Lag Time	Bottle Type
Reagent 1	330 µL	CA1A		120	*
Reagent 2					

* Enter the bottle type: plastic or glass