

LIPASE

Color

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instruments: **EXPRESS 550 / Plus**

Reagent preparation

Reagent 1: use the Reagent A.
 Reagent 2: use the Reagent B.

Instrument settings

Test Name	LIPASE		Test	LIP	
Test Bar Code					
Test Type	Kinetic		Curve Type	Blanked Linear	
Units	U/L		Nr. of Decimal Places	0	
Primary Wavelength	570		Secondary Wavelength		
Read Time/Interval	180		Sample Blank	No	
Factor					
Calibration Interval	168 hours				
Nr. of Calibrators	2		Nr. of Replicates	2	
Low Blank A Limit	-0.100		High Blank A Limit	0.650	
Low A Limit	0.000		High A Limit	2.000	
Low Normal	12		High Normal	48	
Linearity Limit	500		Curve S.D. Limit	8.0	
Sample Volume	4 µL		Test	LIP	
Rerun Dilution Ratio	2		Predilution Ratio	1	
Sample Diluent Bottle Type	*				
Reagent Diluent Bottle Type	*				
	Reagent Volume	Bar Code	Diluent Volume	Lag Time	Bottle Type
Reagent 1	260 µL	LP1A		300	*
Reagent 2	130 µL	LP2A		180	*
			* Enter the bottle type: plastic or glass		