

HEMOGLOBIN A1C

Turbidimetry

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: HITACHI 911 (= HITACHI 904)

Samples

Capillary or venous blood collected by standard procedures and with heparin or EDTA as anticoagulants.
HbA1C in blood is stable 3 days at 15-25°C, 7 days at 2-8°C and 6 months at -20°C.
Freeze once only.

Hemolysate preparation

- The calibrators do not require pretreatment.
1. Bring the reagent A to room temperature.
 2. Pipette into a test tube:

Blood	10 µL
Reagent (A)	1000 µL

3. Shake thoroughly. Avoid the formation of foam. The hemolysate can be used after the solution has changed color from red to brownish-green (approximately 3 minutes).

The hemolysate is stable 4 hours at 15-25°C, 24 hours at 2-8°C and 6 months at -20°C.
Freeze once only.

Reagent preparation

Reagents (A), (B), (C) and (D) are provided ready to use.
HbA1C Standards (S1-S4): Reconstitute with 2.0 mL of distilled water. Stable for 8 hours at 15-25 °C, 2 days at 2-8°C and 3 months at -20°C. Freeze once only.

Hb Reagent 1: Reagent B
HbA1c Reagent 1: Reagent C Reagent 2: Reagent D

Instrument settings

TOTAL HB

TEST	HB	XXX	TEST NAME	HB	UNIT	g/dL
DATA MODE		1: ON BOARD	REPORT NAME	HB		
CONTROL INTERVAL		1000	INSTRUMENT FACTOR (y=aX+b)		a 1.0	b 0.0
EXPECTED VALUE	<S. TYPE 1>		EXPECTED VALUE	<S. TYPE 2>		
AGE	(M)	(F)				
TECHNICAL LIMIT	<S. TYPE 1>			<S. TYPE 2>		
	1.10 - 40					
STD	CONC.	POS.	S.VOL.	PRE.	DIL.VOL.	CODE
(1)	0.0	L	20	0	0	99
(2)	*	...	20	0	0	...
						Blank: Make reagent blank with sodium chloride 154 mmol/L. Calibrator: Standard S4.
TEST	HB		WAVELENGTH (SUB/MAIN)			660 / 570
ASSAY CODE	END POINT - 10 - 0	ASSAY POINT 15 - 0 - 0 - 0	DILUTION			301 / 99
		<S. TYPE 1>				<S. TYPE 2>
S. VOL. (REGULAR)	20 - 0 - 0					
S. VOL. (DECREASE)	20 - 0 - 0					
S. VOL. (INCREASE)	20 - 0 - 0					
ABS. LIMIT	0		0		2: INCREASE	
PROZONE LIMIT	32000				1: UPPER	
REAGENT	R1 230 - 0 - XXX - 28					
	R2 0 - 0 - XXX - 0					
	R3 0 - 0 - XXX - 0					
	R4 0 - 0 - XXX - 0					
CALIB. TYPE	1: LINEAR	2 2 0				
TIME OUT BLANK	0		SD LIMIT			0.1
2 POINT	0		DUPLICATE LIMIT			100
FULL	0		SENSITIVITY LIMIT			0
CHANGE LOT	0		S1 ABS LIMIT	-32000		32000
BOTTLE	0					

HBA1C

TEST	HBA1C		XXX	TEST NAME	HBA1C	UNIT	g/dL	
DATA MODE			1: ON BOARD	REPORT NAME	HBA1C			
CONTROL INTERVAL			1000	INSTRUMENT FACTOR (y=aX+b)		a	1.0	b 0.0
EXPECTED VALUE	<S. TYPE 1>			EXPECTED VALUE	<S. TYPE 2>			
AGE	(M)		(F)					
	* - *		* - *					
TECHNICAL LIMIT	<S. TYPE 1>			<S. TYPE 2>				
	0.05 - 2.50							
STD	CONC.	POS.	S.VOL.	PRE.	DIL.VOL.	CODE		
(1)	0.00	L	10	0	0	999		
(2)	*	...	10	0	0	...		
(3)	*	...	10	0	0	...		
(4)	*	...	10	0	0	...		
(5)	*	...	10	0	0	...		
Blank: Make reagent blank with sodium chloride 154 mmol/L. Calibrators 2-5: Standards S1-S4								
TEST	HBA1C			WAVELENGTH (SUB/MAIN)		700 / 340		
ASSAY CODE	2 POINT END - 10 - 0		ASSAY POINT 15 - 31 - 0 - 0	DILUTION		301 / 99		
	<S. TYPE 1>			<S. TYPE 2>				
S. VOL. (REGULAR)	10 - 0 - 0							
S. VOL. (DECREASE)	10 - 0 - 0							
S. VOL. (INCREASE)	10 - 0 - 0							
ABS. LIMIT	0			0	2: INCREASE			
PROZONE LIMIT	0				2: LOWER			
REAGENT	R1	250 - 0 - XXX - 28						
	R2	0 - 0 - XXX - 0						
	R3	50 - 0 - XXX - 28						
	R4	0 - 0 - XXX - 0						
CALIB. TYPE	3: LOGIT-LOG (4P)		5	4	0			
TIME OUT BLANK			0	SD LIMIT	250			
2 POINT			0	DUPLICATE LIMIT	400			
FULL			0	SENSITIVITY LIMIT	3000			
CHANGE LOT			0	S1 ABS LIMIT	-32000	32000		
BOTTLE			0					

Version 0508

CALCULATION

$$\% \text{HbA1C - IFCC} = \frac{\text{HbA1C (g/dL)}}{\text{Hb (g/dL)}} \times 100$$

The HbA_{1c} percentage in the sample is calculated using the following general formula. The values are traceable to IFCC Reference Method:

The traceable values to Reference Method as described by the US National Glycohemoglobin Standardization Program (NGSP) are calculated using the following general formula:

$$\% \text{HbA1C-NGSP} = 0.915 \times \% \text{HbA1C-IFCC} + 2.15$$