

HEMOGLOBIN A1C

Turbidimetry

The application parameters comprised here constitute a guide to facilitate the validation of our reagents by the instrument. It is advisable to validate the use when there is any change in software or reagent versions.

Instrument: HITACHI 917

Samples

Capillary or venous blood collected by standard procedures and with heparin or EDTA as anticoagulants.

HbA1C in blood is stable 3 days at 15-25°C, 7 days at 2-8°C and 6 months at -20°C. Freeze once only.

Hemolysate preparation

The calibrators do not require pretreatment.

1. Bring the reagent A to room temperature.
2. Pipette into a test tube:

Blood	10 µL
Reagent (A)	1000 µL

3. Shake thoroughly. Avoid the formation of foam. The hemolysate can be used after the solution has changed color from red to brownish-green (approximately 3 minutes).

The hemolysate is stable 4 hours at 15-25°C, 24 hours at 2-8°C and 6 months at -20°C. Freeze once only.

Reagent preparation

Reagents (A), (B), (C) and (D) are provided ready to use.

HbA1C Standards (S1-S4): Reconstitute with 2.0 mL of distilled water. Stable for 8 hours at 15-25 °C, 2 days at 2-8°C and 3 months at -20°C. Freeze once only.

Hb Reagent 1: Reagent B
HbA1c Reagent 1: Reagent C Reagent 2: Reagent D

Instrument settings

HB

ANALIZE		CALIB		
Assay/Time/Point	End Point / 10 / 16 - 0	Calibration Type	Linear	
Wavelength (Second./Main)	660 / 570	Points/Span Point	2 / 2	
Sample Volume (Normal)	20 / 0 / 0	Weight	0	
(Decreased)	10 / 0 / 0	SD Limit	0.1	
(Increased)	20 / 0 / 0	Duplicate Limit % / Abs	10 / 50	
Diluent/Expiration	951 / 99	Sensitivity Limit	-9999 / 9999	
Reagent (R1) T1	230 / 0 / ... / 28	Abs. S1 Limit	-32000 / 32000	
(R2) T2	0 / 0 / ... / 0	OTHERS		
(R3) T3	0 / 0 / ... / 0	<i>Calib. Code</i>	<i>Concentration</i>	<i>Position</i>
(R4) T4	0 / 0 / ... / 0	(1) ...	0.0	...
Abs. Limit	32000 – Increase	(2) ...	*	...
Prozone Limit	-32000 – Lower	<i>Sample volume</i>	<i>Diluted</i>	<i>Diluent</i>
Washing Solution	Detergent 1	(1) 20	0	0
RANGE		(2) 20	0	0
Application N° / Units	... / g/dL	Blank: Make reagent blank with sodium chloride 154 mmol/L.		
Name	Hemoglobin	Calibrator: Standard S4		
Control Interval	1000	... Data entered by the operator		
Instrument Factor	a = 1.0 b = 0.0	(*) Enter S4 Standard value		
Technical Limit	1.10 - 40			
Repetition Limit	1.10 - 40			
References Values Male	... - ...			
Female	... - ...			

HBA1C

ANALYZE		CALIB	
Assay/Time/Point	2 POINT END / 10 / 16 - 34	Calibration Type	Logit-Log (4P)
Wavelength (Second./Main)	700 / 340	Points/Span Point	5 / 4
Sample Volume (Normal)	10 / 0 / 0	Weight	0
(Decreased)	5 / 0 / 0	SD Limit	250
(Increased)	20 / 0 / 0	Duplicate Limit % / Abs	10 / 50
Diluent/Expiration	951 / 99	Sensitivity Limit	-9999 / 9999
Reagent (R1) T1	250 / 0 / ... / 28	Abs. S1 Limit	-32000 / 32000
(R2) T2	0 / 0 / ... / 0		
(R3) T3	50 / 0 / ... / 28	OTHERS	
(R4) T4	0 / 0 / ... / 0	<i>Calib. Code</i>	<i>Concentration</i>
Abs. Limit	32000 – Increase	(1) ...	0.00
Prozone Limit	-32000 – Lower	(2) ...	*
Washing Solution	Detergent 1	(3) ...	*
		(4) ...	*
		(5) ...	*
RANGE			
Application N° / Units	... / g/dL	<i>Sample volume</i>	<i>Diluted</i>
Name	Hemoglobin A1c	(1) 10	0
Control Interval	1000	(2) 10	0
Instrument Factor	a = 1.0 b = 0.0	(3) 10	0
Technical Limit	0.05 – 2.5	(4) 10	0
Repetition Limit	0.05 – 2.5	(5) 10	0
References Values Male	... - ...		
Female	... - ...		
... Data entered by the operator		Blank: Make reagent blank with sodium chloride 154 mmol/L. (*) Calibrators 2-5: Standards S1-S4	

Version 0704

CALCULATION

$$\% \text{HbA1C - IFCC} = \frac{\text{HbA1C (g/dL)}}{\text{Hb (g/dL)}} \times 100$$

The HbA_{1c} percentage in the sample is calculated using the following general formula. The values are traceable to IFCC Reference Method:

The traceable values to Reference Method as described by the US National Glycohemoglobin Standardization Program (NGSP) are calculated using the following general formula:

$$\% \text{HbA1C-NGSP} = 0.915 \times \% \text{HbA1C-IFCC} + 2.15$$