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Reagent 1: Reactivo B.

		A25	A15	
GENERAL	Test name	Hb Total endpoint mon.	Hb Total endpoint mon.	
	Analysis mode	whole blood	whole blood	
	Sample type	g/dL	g/dL	
	Units	increasing	increasing	
	Reaction type	No	No	
	Turbidim. test	1	1	
	Decimals	1	1	
	Replicates	-	-	
	Name of assoc. constituent	-	-	
	PROCEDURE	Type of reading	bichrom.	bichrom.
Volumes	Sample	30	30	
	Reagent 1	170	170	
	Reagent 2	-	-	
	Washing	1.2	1.2	
	Predilution factor	-	-	
	Filters	Main	560	560
		Reference	670	670
	Times	Reading 1	300 s	312 s
		Reading 2	-	-
		Reagent 2	-	-
Postdilution factor		1.2	1.2	
CALIBRATION	Type of calibration	specific	specific	
	Calibrator number	1	1	
	Calibrator replicates	3	3	
	Blank replicates	3	3	
	Calibration curve	-	-	
OPTIONS	Blank with saline solution	yes	yes	
	Blank absorbance limit	-	-	
	Kinetic blank limit	-	-	
	Linearity limit	-	-	

Reagent 1: Reactivo C; Reagent 2: Reactivo D.

		A25	A15	
GENERAL	Test name	HbA1C differential bir.	HbA1C differential bir.	
	Analysis mode	whole blood	whole blood	
	Sample type	g/dL	g/dL	
	Units	increasing	increasing	
	Reaction type	Yes	Yes	
	Turbidim. test	2	2	
	Decimals	1	1	
	Replicates	-	-	
	Name of assoc. constituent	-	-	
	PROCEDURE	Type of reading	monoch.	monoch.
Volumes	Sample	8	8	
	Reagent 1	200	200	
	Reagent 2	40	40	
	Washing	1.2	1.2	
	Predilution factor	-	-	
	Filters	Main	340	340
		Reference	-	-
	Times	Reading 1	180 s	168 s
		Reading 2	480 s	480 s
		Reagent 2	195 s	192 s
Postdilution factor		2	2	
CALIBRATION	Type of calibration	specific	specific	
	Calibrator number	4	4	
	Calibrator replicates	3	3	
	Blank replicates	3	3	
	Calibration curve	decreasing polygonal	decreasing polygonal	
OPTIONS	Blank with saline solution	Yes	Yes	
	Blank absorbance limit	-	-	
	Kinetic blank limit	-	-	
	Linearity limit	-	-	

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1. Karl J et al. Development and standardization of a new immunoturbidimetric HbA_{1c} assay. *Klin Lab* 1993;39:991-996.
2. Tietz NW. *Clinical guide to laboratory tests*, 3rd ed. Saunders Co, 1999.
3. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long term complications in insulin-dependent diabetes mellitus. *N Engl J Med* 1993; 329: 977-986.
4. Jeppson JO, Kodold U, Finke A, Hoelzel W, Hoshino T, Miedema K, Mosca A, Mauri P, Paroni R, Thienpont L, Umamoto M, Weykamp C. Approved IFCC reference method for the measurement of HbA_{1c} in human blood. *Clin Chem Lab Med* 2002;40:78-89.
5. Hoelzel W et al. IFCC reference systems for measurement of hemoglobin A_{1c} in human blood and the national standardization schemes in the United States, Japan and Sweden: a method-comparison study. *Clin Chem* 2004;50:166-174
6. Young DS. *Effects of drugs on clinical laboratory tests*, 4th ed. AACC Press, 1995.
7. Bry L, Chen P, Sacks D. Effects of hemoglobin variants and chemically modified derivatives on assays for glycohemoglobin. *Clin Chem* 2001;41:153-163.
8. Roberts WL, De BK, Brown D, Hanbury MC, Hoyer JD, John WG, Lambert TL, Lundell RB, Rohlfing C, Little RR. Effects of hemoglobin C and S traits on eight glycohemoglobin methods. *Clin Chem* 2002;48:383-385.